

## Sealed IFT Chamber

- **Small sealed chamber suitable for:**
  - liquid-vapor
  - liquid-liquid
  - liquid-solid-vapor
  - liquid-liquid-solid
- **Two internal loops provide temperature stabilization with user supplied circulating bath**
- **Temp. range of below ambient to 200°C depending on bath and fluid**
- **Chamber is sealed with Viton O-rings and rated to 100 PSI**
- **Two 25mm diameter windows provide light passage**
- **Removable sample stage for contact angle work - manual movement**
- **Needle ports at top and bottom for pendant drop and bubble-up work**
- **Internal volume - 22ml**
- **Mounting hardware for FTA1000, FTA200 and FTA188**



IFT Chamber Schematic	IFT Chamber and mounting hardware	IFT Chamber

For more information contact [sales@firsttenangstroms.com](mailto:sales@firsttenangstroms.com)

December, 2009