


December, 2009

FTA1000B Manual USB Drop Shape Analyzer Model B 1A3 310

- > Lowest cost model
- > USB powered
- > Laptop capability
- **>** Drop size, 0.5-40μl
- > Manual drop formation, 2cc micrometer syringe
- Manual drop touch-off, 50mm Z precision movement
- > Large specimen table
- > Camera, 15fps at 640x480 resolution
- > Fixed mag lens, 5mm FOV
- > Fta32 Software for video acquisition and full data analysis
- > High quality mechanical construction



Part No.	
В	B frame and Fta32 software. Site-wide license
B 1	Specimen stage, 125mm sq. platform, 50mm Z adjustment
B 0A	Camera, USB1.1, 1/4" sensor, 15 fps
В 003	Microscope, 100mm WD, fixed mag/fixed focus lens
B 000 3	Backlight, 25mm red LED
B 000 01	Dispenser, 2ml micrometer syringe, 50mm Z movement
В 000 000	Electronic control box not required

Computer sold separately.

For more information contact sales@firsttenangstroms.com